

DINO - LITE *Edge*

AM7025X



Model: **AM7025X** Dino- Lite Edge

SPECIFICATIONS

Model	AM7025X Dino-Lite Edge	Operating System Supported	Window 10, 8, 7, Vista, XP MAC OS 10.9 or later
Interface	USB 2.0	Calibration Function	Yes
Product Resolution	5M pixels (2592 x 1944)	Measurement Function	Yes
Frame Rate	10fps in 5MP / 3MP/ 2MP, MJPEG 25fps in 1.3MP, MJPEG 30fps in VGA, MJPEG	Unit Dimensions	8.0cm (H) x 3.2cm (D)
		Unit Weight	117g

HIGHLIGHTS



High Optical Resolution



5.0 Megapixels



Multi Adaptable Design



Robust Housing

PRODUCT SUMMARY

Adopting with Edge class optics and high-resolution sensor, the 5MP Dino-Eye Edge series AM7025X eyepiece offers remarkable image performance yet provides relatively large field of view than its predecessors. With fine color reproducibility, the AM7025X is capable of capturing finest detail of microscope's image for professional applications.

GALLERY

