

# AM4025X

Catering the needs of digitizing conventional microscope, the Dino-Eye Edge 1.3MP AM4025X provides remarkable image quality for most demanding applications. Adopting with Edge-class optics and advanced image sensor, the AM4025X delivers sharp images with authentic color yet to provide large FOV for revealing more details.



Overview	Specification	Gallery
----------	---------------	---------

Model	AM4025X Dino-Eye Edge
Interface	USB 2.0
Product Resolution	1.3M pixels (1280x1024)
Frame Rate	15fps in 1.3MP, MJPEG 30fps in VGA, MJPEG
Operating System Supported	Windows 10, 8, 7, Vista, XP Mac OS 10.12 ~ 10.14
Calibration Function	Yes
Measurement Function	Yes
Unit Dimension	9.2cm (H) x 3.2cm (D)
Unit Weight	80g